

## Raman

Software for Raman spectra analysis.

### Open Single Spectrum Files

Opens one or multiple spectrum files.

Expected format of file: Each line has two values separated by space. First value is x coordinate, second value is y coordinate.

### Open Multi Spectrum Files

Opens one or multiple spectrum files.

Expected format of file:

First line contains x coordinates (separated by space).

Second line contains y coordinates (separated by space) of first spectrum.

Third line contains y coordinates (separated by space) of second spectrum.

### Zoom

User can Zoom In and Zoom Out just scrolling mouse middle button. Dragging mouse middle button moves the view.

### Zoom Windows

Draw rectangle to see specific part of spectrum.

### Zoom To Original Size

Shows whole spectrum.

### Baseline Correction

Click at multiple positions (at least 4) on spectrum to create baseline. You can remove any point by clicking right mouse button and selecting context menu *Remove Closest Point*.

### Peak Analysis

Select border points on baseline which represents the peak. You can remove any point by clicking right mouse button and selecting context menu *Remove Closest Point*.

### Cut Off

Enables to cut off left and right part of spectrum.

### Format of file with exported peaks

Spectrum\_Name,Peak\_Number,Peak\_Start,Peak\_End,Peak\_Height,Peak\_Position,Area