

Raman

Program for Raman spectra analysis.

Baseline Correction

Select on spectrum baseline points. Select at least 4 points to see baseline. You can remove any point by clicking right mouse button and selecting in context menu to remove the point.

Peak Analysis

Select border points on baseline which represents the peak. You can remove any point by clicking right mouse button and selecting in context menu to remove the point.

Zoom

User can Zoom In and Zoom Out just scrolling mouse middle button. Dragging mouse middle button moves the view.

Format of file with exported peaks

File has following columns separated by tabulator:

Peak number

Spectrum name

Peak height

x coordinates of left peak border

x coordinates of right peak border

x coordinates of peak